OGI-IT LIMS

a visual introduction

What is a LIMS

LIMS - **L**aboratory **I**nformation **M**anagement **S**ystem

a **C**omputer **A**pplication to **m**anage a **L**aboratory

Information is represented on a **C**omputer as a **d**ata **m**odel

the **d**ata is transformed (sorted, queried, aggregated) so that
management decisions can be made

Data Model

The **data model** is tightly coupled to the „place“ where the data is stored, the **database**.

- **relational** database - many **tables**, each table can have many **columns**, each column can hold only one **type** of data, each cell holds one **value**
- **document oriented** database - many **collections**, each collection consists of many similar **documents**, each document has a **complex structure** (objects, attributes, array's etc.)

Data Model of a Laboratory

Some entities in a Laboratory are:

- Studies (Studien)
- Samples (Studienproben)
- Operating Procedures (Laborvorschriften)
- Batches (Laborarbeiten)
- Equipment (Laborgeräte)
- External Materials (externe Materialien) i.e. chemicals, pure substances, matrices, etc.
- etc.

in a **document oriented** LIMS these are **collections**

Collections and Documents

The screenshot displays the OGI-IT LIMS v.0.0.4 interface. On the left, a sidebar shows 'Verfügbare Sammlungen' (Available Collections) with a tree view including 'Studienproben', 'Laborvorschriften', 'Laboraufträge', 'Laborarbeiten', 'Laborgeräte', 'externe Materialien', and 'Laborprodukte'. The main area shows a list of external materials under 'externe Materialien' with columns for 'name', 'typ', and 'Bezeichnung'. A detailed view of a document for 'H2O' is shown, including fields for 'UUID', 'Name', 'Typ', 'Validität Temperatur', 'Validität Verfall', 'Berechnung Werte', and 'Daten'. The 'Daten' section includes 'Hersteller', 'Katalog-/Bestellnummer', 'Batch-/Charge-/Serien-No', and 'CAS-Registrierungsnummer'. A 'Substance Details' panel on the right shows 'CAS Registry Number: 7732-18-5' and 'CA Index Name: Water', along with a list of synonyms.

- the LIMS consists of „collections“
- „collections“ contain „documents“
- „documents“ are arbitrary complex structures of attributes like „strings“, „numbers“ but also „array’s“ of data or „links“

Customization

~~Laboratory~~ **Information** **Managent** **System**

IMS the foundation of a **LIMS**

LIMS - a customized IMS

individually define the **Laboratory**

Business Process

business process - a **sequence** of related, **structured activities** or **tasks** that produce a specific service or product (http://en.wikipedia.org/wiki/Business_process)

LIMS workflow - **all** possible **states** of a **document** in the **data model** of the **Lab**

signature - the signee „**guarantees**“ the **content** of a **document** at the time of signing

Audit Trail - the **guaranteed** current and previous **states** of a **document**

Audit Trail - the **array** of „**signatures**“ in a **document**

LIMS document workflow - the **array** of „**signatures**“ in a **document**

LIMS deleted documents - collection to store any kind of removed document

Name (1) datetest2
Absender (1) xxxx
Lieferdatum erwartet (1) 7.5.2015, 00:00:00
Frächter (2) DHL
Frachtbrief Nr. (2) Tracking Nummer z.B. HWB-Nummer.

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Unterschriften Mappe

Status der Unterschriften received

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
properties of derived probes are calculated automatically

Analyst Companion App

The OGI-IT LIMS Companion App is started inside of the SCIEX Analyst Software, a widely used LC/MS/MS (Liquid chromatography–mass spectrometry/ mass spectrometry) instrument control software. It allows both import of „Acquisition Batches“ from the LIMS and export of „Quantitation Results“ back into the LIMS.

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

Standard-Drucker

Etiketten drucken

support for template network label printers

Visual Feedback

- image-based barcode reading
- get visual feedback when scanning

OGI-IT LIMS v0.6.3

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- IB2

Studie	Name	Status	Tags	Lagerort
xxx	IA1	OK	Timo1	TK 26
				TK 26
				TK 26
				TK 22
				TK 26
				TK 26
				TK 26
				TK 26
				TK 27
				TK 26
				TK 26
				TK 22

Kontext

Benutzer

Standard-Scanner2
Standard-Scanner3

50036656
FC13780848
FC13780848
FC13780836
FC13780812
FC13780823
FC13780836
FC13780848
FC13780848
FC13780824
FC13780836
50036656
50036656

Löschen

Laborprodukte: QK1

Name IB2
UID 21ada87c-ff4d-43be-b851-39493e770e42
Studie xxx
Status OK

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
OGI-IT LIMS Reporting Worker v.0.5.7
Probenverzeichnis

- template based Reporting
- PDF
 - CSV
 - JSON

Quality

General

- **ISO-9000 ff** - http://en.wikipedia.org/wiki/ISO_9000 family of quality management systems standards is designed to help organizations ensure that they meet the needs of customers and other stakeholders while meeting statutory and regulatory requirements related to a product

Specific to Laboratory

designed with this in mind

- **ISO-15189** - http://en.wikipedia.org/wiki/ISO_15189 based on ISO/IEC 17025 (General requirements for the competence of testing and calibration laboratories) and ISO 9001, it is a unique document that takes into consideration the specific requirements of the medical environment and the importance of the medical laboratory to patient care
- **GAMP 5** - http://en.wikipedia.org/wiki/Good_Automated_Manufacturing_Practice a set of guidelines for manufacturers and users of automated systems in the pharmaceutical industry (ISPE - International - EU)
- **CFR Part 11**- http://en.wikipedia.org/wiki/Title_21_CFR_Part_11 defines the criteria under which electronic records and electronic signatures are considered to be trustworthy, reliable and equivalent to paper records (FDA -US)

Thank you!